



OXYGEN CONCENTRATOR C5A NI

ISO
9001
Certified

ISO
13485
Certified

You know we care..



OXYGEN CONCENTRATOR

C5A N1

SPECIFICATIONS :

- **Power supply:** AC230V, 50Hz/AC220V, 60Hz/AC120V, 650Hz; Power: 390VA
- **Sound level:** $\leq 46\text{dB (A)}$
- **Maximum recommended flow:** 5L/min
- **Flow Range at Outlet Pressure:** 0.5~5L/min
- **Oxygen Concentration:** 0.5~5L/min, 93% \pm 3%
- **Output Pressure :** 38kPa \pm 5kPa;
- **Release Pressure by machine operation:** 250kPa \pm 50kPa;
- **Net Weight:** 16.0kg, Gross Weight: 18.5kg;
- **Product Dimension:** 33 x 26 x 54 cm
- **Package Dimension:** 37 x 43 x 70 cm

SAFETY SYSTEM :

- Current overload or line surge shutdown.
- High temperature compressor shutdown.
- High pressure alarm shutdown.
- Low pressure alarm shutdown.
- Low Oxygen Concentration alarm.



OUR SOLUTIONS

- Pneumatic Waste & Laundry Collection System (PWLCS)
- Bio-Medical Waste Treatment System (ISS)
- Pneumatic Tube System (PTS)
- Central Sterile Supply Department (CSSD)
- Modular OT / ICU / Emergency Room / OT Integration
- Medical Gas Pipeline System (MGPS)
- Oxygen Generation System (OGS)
- Nurse Call System (NCS)
- Automated Inventory Management Systems (AIMS)
- Burns Ward & Equipment
- Hospital Furniture / OT Table and OT Lights / OT Equipment
- Curtain Track / IV Track / Wall Guard
- Smart Bin System (SBS)
- Electric Track Vehicle (ETV)
- Automated Pharmacy Dispensing System (APDS)
- Hospital Consultancy & Designing on Turnkey Basis

Imported & Marketed by



UNG Tower, 3rd Floor, Building No. 2, Central Market,
Punjabi Bagh West, New Delhi - 110026, INDIA

Phones : (+91-11) 4545 9999 (Hunting Line)

E-mail : sales@medfreshe.com, info@medfreshe.com

Website : www.medfreshe.com

 MedFreshe |  MedFreshe

 MedFreshe |  unisonnarulagroup

Narula's are proud to be members of

